



Amplivox Screening Audiometer Comparison Matrix

User Features	Amplivox 116	Amplivox 170
Portable	✓	✓
PC connectivity (USB)	✓	✓
Thermal printer output	✓	✓
Manual test	✓	✓
Automatic test	-	✓
Internal memory	✓	✓
Communication channel to patient	✓	✓
Patient response switch	✓	✓
Mains power adapters	✓	✓
Technical Features		
Frequency range	125-8K	125-K
Attenuation range	-10 to 100 max	-10 to 100 max
Attenuation step size (dBHL)	5	5
Pure tone signal	✓	✓
Warble tone signal	✓	-
Continuous signal	✓	-
Pulsed signal	✓	-
Dimensions	10.62" x 6.5" x 2.36"	10.62" x 6.5" x 2.36"
Weight	1.527 lbs	1.565 lbs
External PSU	✓	✓
Optional battery power	✓	-